

## Methenyl-THF Cyclohydrolase Assay

Assay monitors the hydrolysis of 5,10-methenyl-THF to 10-formyl-THF by following the decrease in absorbance at 355 nm. The substrate, 5,10-methenyl-THF, has a high non-enzymatic buffer-catalyzed hydrolysis rate, which must be subtracted from the experimentally-measured rate. This non-enzymatic hydrolysis can be minimized by using potassium maleate as the buffer, which has a 40-fold lower rate than phosphate buffer (Rabinowitz, J. C. (1960) in *The Enzymes* (Boyer, P. D., ed) Vol. 2, 2nd Ed., pp. 185-252, Academic Press, New York).

### Procedure:

Prepare fresh cocktail:

5 ml 1 M K•maleate (pH 7.4) (use KOH to adjust pH of maleic acid)

34  $\mu$ l 2-mercaptoethanol (14.6 M stock)

14.2 ml H<sub>2</sub>O

19.375 ml = ~50 assays at 375  $\mu$ l/tube

Leaves room for 100  $\mu$ l sample/buffer + 25  $\mu$ l substrate = 0.5 ml total reaction volume. Samples should be diluted in the buffer cocktail to minimize non-enzymatic hydrolysis.

Substrate: stock solution is ~2 mM. Dilute 10  $\mu$ l into 1 ml 0.36 N HCl and measure A<sub>350</sub> to calculate actual concentration:

$$A_{350}/24.9 \text{ mM}^{-1} \times 100 = \text{mM}$$

Final assay concentrations are:

0.2 M K•maleate (pH 7.4)

20 mM 2-mercaptoethanol

~100  $\mu$ M (R,S)5,10-methenyl-THF

Blank the spec at 355 nm with 20  $\mu$ l substrate in 0.5 ml 20 mM 2-mercaptoethanol/0.1 N HCl (13.5  $\mu$ l BME into 10 ml 0.1 N HCl)

### Assay:

- add 375  $\mu$ l cocktail to cuvette
- add 100  $\mu$ l sample/buffer
- initiate rxn with 25  $\mu$ l substrate stock; cap cuvette with patch of parafilm and invert 3X to mix.
- place in spec and start recording

Data reduction: calculate  $\Delta A_{355}$  over 15-60 sec (may need to ignore first 10-15 sec to get stable slope), then convert to  $\Delta A_{355}/\text{min}$  by multiplication.

Need to do a no-enzyme blank to determine the non-enzymatic rate of substrate hydrolysis – this will be subtracted from each enzyme sample. Typical rate for no-enzyme blank is 0.03-0.06  $\Delta A_{355}/\text{min}$

For the SpectraMax, find the Cyclo file in the Protocol Folder of the SoftMax folder. Choose CuvetteSet; Kinetics; Vmax. To blank the instrument, read the cuvette with substrate + BME/HCl using the REF button. Samples can then be read using the READ button. They can each be named using the Template function, and a table can be printed out with all the data. The software reports the the V<sub>max</sub> as change in milli A<sub>355</sub>/min [i.e.  $\Delta(A_{355} \times 10^{-3})$  per min].

$$\frac{1 \mu\text{M}}{0.0249 A_{355}} \square \square A_{355}/\text{min} = \frac{1 \mu\text{mol}}{l} = \frac{\text{nmol}}{\text{ml}} \square 0.5 \text{ml} = \text{nmol product}/\text{min}$$

$$\boxed{\Delta A_{355}/\text{min} \times 20.1 = \text{nmol}/\text{min} = \text{mU}}$$