

## 10-Formyl-THF Synthetase Assay

Assay monitors the ATP-dependent formylation of THF to 10-formyl-THF. 10-formyl-THF is converted non-enzymatically to 5,10-methenyl-THF by the addition of acid to stop the reaction. The absorbance at 350 nm is used to calculate the amount of product formed, using  $\epsilon_{350} = 24900 \text{ M}^{-1}\text{cm}^{-1}$

### Stocks

1 M Tris-SO<sub>4</sub>, pH 7.5

1 M KCl

1 M NH<sub>4</sub>formate, pH 7.5

[titrate 1 M formic acid (prepared from 88% formic acid = 23.4 M) with NH<sub>4</sub>OH to pH 7.5]

0.1 M MgSO<sub>4</sub>

50 mM K•ATP, pH 7.5

(0.138 g/5 ml; adjusted to pH 7.5 with KOH)

10 mM THF (prepared as described)

Final assay concentrations are:

25 mM Tris-SO<sub>4</sub>/100 mM KCl/10 mM MgSO<sub>4</sub>/100 mM NH<sub>4</sub>formate/5 mM K•ATP/2 mM THF

### Procedure:

Prepare assay cocktail (can be stored at 4°C for up to 1 month):

1.25 ml Tris-SO<sub>4</sub>

5 ml MgSO<sub>4</sub>

5 ml KCl

5 ml NH<sub>4</sub>formate

5 ml ATP

13.75 ml H<sub>2</sub>O

35 ml = 100 assays at 350 µl/tube

Add 350 µl cocktail + 100 µl 10 mM THF to each tube and preincubate at 37°C for 1 min.

Initiate reaction with 50 µl enzyme/buffer = 0.5 ml total reaction volume. Incubate 10 min at 37°C.

Stop reactions with 1 ml 0.36 N HCl. Vortex, let stand at RT 5 min.

Read absorbance at 350 nm. Include a no enzyme blank. Read blank against water first – if THF is good, should read < 0.05. Re-blank spec against this, and read samples.

Data reduction:

$$A_{350} \times \frac{1 \mu\text{M}}{0.0249 A_{350}} \times \frac{1}{10 \text{ min}} = \frac{1 \text{ } \mu\text{mol}}{\text{l} \cdot \text{min}} = \frac{\text{nmol}}{\text{ml} \cdot \text{min}} \times 1.5 \text{ ml} = \text{nmol product / min} = \text{mU}$$

$$\boxed{A_{350} \times 6.02 = \text{nmol/min} = \text{mU}}$$